

## User Guide: Basement Membrane Annotation & GBM Profiling Macro (ImageJ)

### Overview

This macro allows users to annotate basement membrane boundaries in .tif images using the freehand tool. It then:

- Saves annotations as ROI zip files,
  - Creates a binary mask,
  - Calculates basement membrane thickness (via distance map),
  - Skeletonizes the membrane,
  - Computes width values,
  - Saves the results as heatmaps and CSV profiles.
- 

### Requirements

- ImageJ or Fiji
  - Images in .tif format
  - Calibration set in microns (optional, but recommended)
- 

### Input Parameters (Prompted Automatically)

When the macro runs, it will prompt the user to specify:

1. **Input Directory** – Folder containing .tif images to process.
  2. **Output Directory** – Folder where results will be saved.
  3. **File Suffix** – Usually .tif (default).
- 

### Step-by-Step Instructions

1. **Launch ImageJ/Fiji** and open the macro (Plugins > Macros > Run...).
2. **Select the input/output directories** and confirm the file suffix.
3. **For each image:**
  - You'll be prompted to annotate the **basement membrane** using the **Freehand tool**.
  - After drawing, click **OK**.
  - You'll be asked if you want to annotate **more regions in the same image**.
    - Click **Yes** to draw more, or **No** to finish the image.

- Annotations will be saved as [filename]\_BM.zip.
4. The macro will then:
- Combine all ROIs
  - Generate a binary mask and compute the **distance map (heatmap)**.
  - Skeletonize the membrane.
  - Calculate **GBM width** as pixel/micron values.
  - Save:
    - Heatmap image: [filename]\_HeatMap.tif
    - Final GBM width image: BM Width Image
    - Width profile as CSV: [filename]\_GBMprofile.csv
5. You will be asked if you want to proceed to the **next image**.
- 

## Output Files

For each input image, the following files are saved:

- \*\_BM.zip → Basement membrane ROI annotations
  - \*\_HeatMap.tif → Distance map of annotated membrane
  - \*\_GBMprofile.csv → Width values of basement membrane
  - BM Width Image (in memory only unless saved manually)
- 

## Notes

- Ensure that the **pixel size is set correctly** in the image metadata; otherwise, the macro defaults to pixels.
- For precise GBM thickness (in microns), both pixel width and height must match and unit must be "microns".
- The macro performs basic file management (closes images after processing).
- Skeletonized membranes are scaled to generate a pixelwise width profile.